

## Pengaruh asam askorbat, karbon aktif, dan durasi pencahayaan terhadap pencokelatan eksplan dan regenerasi tunas kultur in vitro pisang Barangan (*Musa acuminata* L.)

Kusuma Kariyana, author

Deskripsi Lengkap: <https://lib.ui.ac.id/detail?id=20353549&lokasi=lokal>

---

### Abstrak

The research was conducted to determine the effect of ascorbic acid (50 mg/l, 100 mg/l, 200 mg/l) and activated charcoal (0.5 g/l, 1 g/l, 2 g/l) independently with different light duration (4 weeks in darkness, 2 weeks in darkness followed by 2 weeks in 16 hours light and 4 weeks in 16 hours light) on shoot regeneration. Explants of banana cultivar Barangan (*Musa acuminata* L.) were planted on MS basal media supplemented with 1.6 mg/l IAA, 4.0 mg/l BAP and cultured for 4 weeks. After 4 weeks, explant browning level was evaluated. Explants were then cut vertically into two pieces and planted on the same media to induce shoot regeneration. After 4 weeks in shoot regeneration media, number of shoot, colour of shoot and height of shoot were evaluated. MS media supplemented with 1.6 mg/l IAA and 4.0 mg/l BAP without ascorbic acid and activated charcoal in darkness for 4 weeks was the most suitable media for shoot regeneration. The shoot regeneration gave average of 10,4 shoots per explant.

.....

The research was conducted to determine the effect of ascorbic acid (50 mg/l, 100 mg/l, 200 mg/l) and activated charcoal (0.5 g/l, 1 g/l, 2 g/l) independently with different light duration (4 weeks in darkness, 2 weeks in darkness followed by 2 weeks in 16 hours light and 4 weeks in 16 hours light) on shoot regeneration. Explants of banana cultivar Barangan (*Musa acuminata* L.) were planted on MS basal media supplemented with 1.6 mg/l IAA, 4.0 mg/l BAP and cultured for 4 weeks. After 4 weeks, explant browning level was evaluated. Explants were then cut vertically into two pieces and planted on the same media to induce shoot regeneration. After 4 weeks in shoot regeneration media, number of shoot, colour of shoot and height of shoot were evaluated. MS media supplemented with 1.6 mg/l IAA and 4.0 mg/l BAP without ascorbic acid and activated charcoal in darkness for 4 weeks was the most suitable media for shoot regeneration. The shoot regeneration gave average of 10,4 shoots per explant.